



**Agilent Technologies**

Dear Scientist,

Thank you for downloading the Agilent 5977B GC MS Application Notes.

As part of Agilent Technologies and Gulf Bio Analytical group of companies commitment to helping you select the best methods for your applications, we are providing this handy 22 page document which cover the following applications :

- Pesticide Analysis
- Volatiles Analysis
- Semivolatiles Analysis
- Drug Analysis

## **Gulf Bio Analytical Group of Companies**

*Your trusted partner for all analytical needs in UAE, KSA, Oman, Qatar and Bahrain*


Established in 1999 in Dubai, United Arab Emirates, Gulf Bio Analytical LLC (GBA) is now a recognized leader in providing comprehensive chromatography and mass spectrometry laboratory solutions from Agilent Technologies to chemistry and life sciences research laboratories in the Middle East, India and S.E. Asia. From our headquarters in Dubai and all our regional offices we are able to efficiently serve the growing laboratory clientele in this region.

We hope that these application notes prove to be a useful resource enabling you to meet your analytical goals and objectives.

Thanking you,

Gulf Bio Analytical Group  
[contactus@gulfbioanalytical.com](mailto:contactus@gulfbioanalytical.com)





# Identify More Pesticides Faster Using the Agilent High-Efficiency Source

The Agilent Pesticide DRS Screening GC/MSD Analyzer

## Application Note

Food Testing & Agriculture

### Authors

Melissa Churley and Bruce Quimby  
Agilent Technologies, Inc.

The Agilent Pesticide DRS Screening GC/MSD Analyzer, based on the Agilent 7890 GC and the Agilent 5977B GC/MSD, delivers fast screening and quantitation of large numbers of pesticides and endocrine disruptors in a single run. Deconvolution reporting software and a retention-time-locked database of pesticides and endocrine disruptors accelerates reporting and increases the number of targets screened. When configured with the 5977B GC/MSD and a high-efficiency source (HES), the analyzer identifies a greater number of pesticides while reducing analysis time.



**Agilent Technologies**

## Identify More Pesticides with the High-Efficiency Source

Routine analysis of residues in environmental and food samples requires confident identification and low-level detection along with fast reporting times. The GC/MSD pesticide analyzer addresses each of these needs through Deconvolution Reporting Software (DRS), which uses the NIST AMDIS program [1]. The Pesticides and Endocrine Disruptors database [2] accelerates reporting time, and increases the number of targets screened. In addition, proprietary capillary flow technology column backflush shortens cycle time, reduces chemical background, and optimizes uptime.

The 5977B GC/MSD and HES improve screening capability by increasing the number of ions created in the source and transferred into the quadrupole analyzer. More ions deliver more signal and, thus, better sensitivity. This increase in response translates into more targets found during the screening process with good library matches. Positive identification in food samples at 10 ng/g detection levels is now possible using the full-scan mode.

We demonstrated this in a DRS analysis of tomato extract that was spiked with over 200 pesticides at concentrations of 10 and 100 ng/g. This equates to injection of 10 and 100 pg of each pesticide, respectively. At the 10 ng/g level, 38 target compounds were identified using the HES as compared to none using the extractor source. Almost twice as many targets were identified at the 100 ng/g level (Table 1). Figure 1 shows an example AMDIS analysis for the target flusilazole, including raw and extracted spectra with library matching for the component.

## Conclusions

The Agilent Pesticide DRS Screening GC/MSD Analyzer delivers faster and more accurate screening of pesticides when configured with the Agilent 5977B GC/MSD and HES. When combined with Deconvolution Reporting Software, positive identification in full-scan mode for many targets in food at a concentration of 10 ng/g is possible.

Table 1. Number of AMDIS targets identified in tomato spiked at 10 and 100 ng/g using the extractor source (EXR) and HES (MMF = 80). The amount of pesticide injected was 10 and 100 pg, respectively. At the 10 ng/g level, 38 target compounds were identified using the HES as compared to zero using the extractor source. Almost twice as many targets were identified at the 100 ng/g level (Table 1). The NIST hit number breakdown (distribution) is given for categories 1st, 2nd, and  $\geq 3$  hit. Identified targets that were not spiked into tomato but had an AMDIS match score  $\geq 80$  and NIST hit no.  $\leq 3$  are also listed. Tuning conditions for each source are in parentheses.

	EXR (atune)		HES (autotune)	
	10 ng/g	100 ng/g	10 ng/g	100 ng/g
Number of targets with AMDIS score $\geq 80$	0	91	38	164
<b>Distribution of NIST hits</b>				
1st hit	0	63	26	144
2nd hit	0	12	7	14
$\geq 3$ rd hit	0	16	5	6
Not spiked, $\leq 3$ rd hit	2*	4**	2*	8***

\* Diethyl phthalate and benzophenone

\*\* Benzilamide, benzophenone, quitozene metabolite (pentachlorophenyl methyl sulfide), indoxacarb, and dioxacarb decomposition product [phenol, 2-(1,3-dioxolan-2-yl)-]

\*\*\* Diethyl phthalate, benzophenone, fonofos, phenol, phthalic acid, di(oct-3-yl) ester, phthalimide, quitozene metabolite (pentachlorophenyl methyl sulfide), indoxacarb, and dioxacarb decomposition product [phenol, 2-(1,3-dioxolan-2-yl)-]




Figure 1. Analysis of 10 pg flusilazole in tomato using AMDIS. A) Overlay of extracted ion  $m/z$  233 (red) and TIC (black); B) raw spectrum; C) extracted spectrum for the component; D) library spectrum, AMDIS match factor = 84. The reported NIST reverse match score is 73.

## Acknowledgement

The authors wish to thank Nathan Contino.

## References

1. Anon. NIST Standard Reference Database 1A, NIST/EPA/NIH Mass Spectral Library (NIST 14) and NIST Mass Spectral Search Program (Version 2.2), User's Guide. National Institute of Standards and Technology, U.S. Department of Commerce, Gaithersburg, MD, USA.  
<http://www.nist.gov/srd/upload/NIST1aVer22Man.pdf>
2. Wylie, P. L. *Screening for 926 Pesticides and Endocrine Disruptors by GC/MS with Deconvolution Reporting Software and a New Pesticide Library*; Application note, Agilent Technologies, Inc. Publication number 5989-5076EN, **2006**.

## For More Information


These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015  
Printed in the USA  
October 19, 2015  
5991-6355EN



# Improved Volatiles Analysis Using Static Headspace, the Agilent 5977B GC/MSD, and a High-efficiency Source

## Application Note

Environmental

### Authors

Peter Gautschi and Harry Prest  
Senior Application Scientist and  
Senior R&D Scientist  
Agilent Technologies, Inc.

### Introduction

A preliminary assessment of environmental volatiles in water using static headspace analysis was performed using the Agilent 7697A Headspace Sampler, Agilent 7890B GC, and the Agilent 5977B GC/MSD system with high efficiency source (HES). The revolutionary design of the HES produces a higher ion current yield for many compounds (greater sensitivity), which allows flexible approaches to sample analysis such as lowering detection limits, reducing sample size, speeding up analysis, and so forth.

Volatile organic analysis (VOA) by static headspace is a widely applied technique, but it is also particularly challenging as the response factors for the many potential analytes vary widely. This application note presents a survey of select compounds of environmental interest as an indication of what may be achieved with the 5977B GC/MSD in this approach.

Analysis was performed in selected ion monitoring mode of a mixture of VOA compounds spiked into reverse osmosis (RO) water over a calibration range of 0.02–20 µg/L. Replicate injections were made at 0.04 µg/L to assess the method detection limits (MDL). A study of replicates of local tap water was used to demonstrate long-term stability for some naturally occurring compounds.



**Agilent Technologies**

## Materials and Methods

### Sample preparation

A 10 mL volume of water was added to each vial. Aliquots of Stock Standard (prepared in methanol) and Stock Internal Standard (prepared in methanol) were spiked into the solution and the vial sealed. Standards were prepared at 0.02, 0.05, 0.1, 0.2, 0.5, 1, 10, and 20 µg/L.

Table 1 presents a summary of the instrumental acquisition conditions.

Table 1. Headspace Conditions and MS Conditions

Headspace parameters	Agilent 7697A Headspace Sampler	GC parameters	Agilent 7890B GC
<b>Instrument settings</b>		<b>Inlet</b>	
Loop size	1 mL	Inlet type	Split/Splitless Inlet (SSL)
Transfer line type	Fused Silica, deactivated (p/n160-2535-5)	Mode	Split
Transfer line diameter	0.53 mm	Inlet liner	Straight, 2 mm id, 250 µL (p/n 5181-8818)
HSS-GC coupling	Transfer line interface (G3520A)	Heater	125 °C
Carrier control	GC Instrument	Column flow	1.5 mL/min constant flow
Pressurization gas	Helium	Total flow	25 mL/min
Vial standby flow	20 mL/min	Septum purge flow	1.0 mL/min
<b>Temperature settings</b>		Gas saver	OFF
Oven temperature	75 °C	Split ratio	15:1
Loop temperature	75 °C	Split flow	22.5 mL/min
Transfer line temperature	110 °C	<b>Oven</b>	
Transfer line interface (Aux1)	115 °C	Column	Agilent VF-624 MS
<b>Timing settings</b>		Column dimensions	60 m × 0.25 mm, 1.4 µm
Vial equilibration time	12 minutes	Equilibration time	0.25 minutes
Injection duration	0.3 minutes	Temperature program	32 °C (2 minutes), 12 °C/min to 220 °C (5 minutes)
GC cycle time	30 minutes	<b>Mass Selective Detector parameters Agilent 5977B</b>	
<b>Vial and loop settings</b>		Type	High Efficiency Source (HES EI)
Vial size	20 mL	Source temperature	300 °C
Vial shaking	Level 7	Quad temperature	150 °C
Fill pressure	10 psi	Transfer line temperature	280 °C
Fill time	0.2 minutes	Tune file	HES Auto Tune (HES_Atune.u)
Loop ramp rate	20 psi/min	Acquisition type	SIM (see Table 2)
Loop final pressure	7 psi	Solvent delay	3.95 minutes
Loop equilibration time	0.01 minutes	Gain factor	3
Post injection purge	100 mL/min for 2 minutes		
Leak check	Default, 0.2 mL/min		
Mode	Single Extraction		

## Results and Discussion

Table 2 gives the results of an MDL study performed at 0.04 µg/L, with nine replicate analyses. Note that all MDLs are below 0.025 µg/L or 25 ppt with the exception of two compounds, which have MDLs below 30 ppt. The majority of compounds produce MDLs below 0.015 µg/L, including some compounds with relatively low response.

Table 2. A Summary of Selected Compounds Acquired Stating Retention Times in Minutes, Target Quantitation Ion (SIM), and MDL Calculated for Nine Replicates at 0.04 µg/L

Name	RT	Quant ion	MDL	Name	RT	Quant ion	MDL
Vinyl chloride	4.934	62	0.004	1,2-Dibromoethane	13.427	106.9	0.006
Bromomethane	5.611	93.9	0.003	Chlorobenzene	13.969	112	0.015
Chloroethane	5.806	64	0.003	Ethylbenzene	14.03	91	0.014
1,1-Dichloroethene	7.007	95.9	0.008	1,1,1,2-Tetrachloroethane	14.049	130.9	0.005
<i>trans</i> -1,2-Dichloroethene	8.007	95.9	0.009	<i>o</i> -Xylene	14.664	91	0.018
1,1-Dichloroethane	8.554	63	0.004	Styrene	14.683	104	0.015
<i>cis</i> -1,2-Dichloroethene	9.19	95.9	0.011	Bromoform	14.975	170.8	0.006
2,2-Dichloropropane	9.208	77	0.013	1,1,2,2-Tetrachloroethane	15.45	82.9	0.041
Bromochloromethane	9.47	127.8	0.004	1,2,3-Trichloropropane	15.567	110	0.007
1,1,1-Trichloroethane	9.769	96.9	0.005	Bromobenzene	15.573	155.9	0.017
1,1-Dichloro-1-propene	9.921	75	0.012	<i>n</i> -Propylbenzene	15.63	91	0.017
Carbon tetrachloride	9.94	116.9	0.003	2-Chlorotoluene	15.768	91	0.016
Benzene * (blank issue)	10.165	78	0.009	1,3,5-Trimethylbenzene	15.84	105	0.018
1,2-Dichloroethane	10.202	62	0.006	4-Chlorotoluene	15.914	91	0.018
Trichloroethene	10.848	129.9	0.009	<i>tert</i> -Butylbenzene	16.225	134	0.017
1,2-Dichloropropane	11.165	63	0.005	<i>sec</i> -Butylbenzene	16.499	105	0.016
Dibromomethane	11.275	173.8	0.006	4-Isopropyltoluene	16.67	119	0.017
Bromodichloromethane	11.421	82.9	0.005	1,3-Dichlorobenzene	16.719	145.9	0.020
<i>cis</i> -1,3-Dichloropropene	11.89	75	0.014	1,4-Dichlorobenzene	16.841	145.9	0.023
<i>trans</i> -1,3-Dichloropropene	12.506	75	0.013	<i>n</i> -Butylbenzene	17.194	134	0.020
1,1,2-Trichloroethane	12.762	96.9	0.011	1,2-Dichlorobenzene	17.316	145.9	0.021
Tetrachloroethene	12.884	163.8	0.009	1,2-Dibromo-3-chloropropane	18.334	154.9	0.010
1,3-Dichloropropane	12.963	76	0.009	1,2,4-Trichlorobenzene	19.493	179.9	0.028
Dibromochloromethane	13.238	126.8	0.004	Hexachlorobutadiene	19.651	224.8	0.006

\*Blanks showed some low-level contamination for benzene.



Figure 1 shows an example of the linearity achieved over the concentration range 0.02 to 20 µg/L for several representative compounds. Figures 2 and 3 show the system stability for a few compounds of interest. It is noteworthy that these preliminary results do not benefit from the use of internal standard corrections (that is, external standard calibration, and so forth), which can be expected to greatly improve all aspects of the analysis.

## Conclusions

These preliminary results suggest a significant improvement in detection limits is possible in VOA applications through the HES of the Agilent 5977B GC/MSD. The signal improvement provided is not complicated by interferences, and results in clear enhancements in detection.




Figure 1. Linearity to 20 µg/L (external standard) for *trans*-1,2-dichloroethene (A), 1,2-dibromoethane (B), and bromoform (C).




Figure 2. Overlay of the EIC for nine replicate injections of vinyl chloride at 0.04 µg/L.




Figure 3. Response as peak areas for replicate injections (n = 20) of selected compounds showing system stability in time: benzene ~0.01 ppb, dibromomethane 1.5 ppb, and bromochloromethane 0.8 ppb.

## **For More Information**

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)


Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2016  
Printed in the USA  
January 7, 2016  
5991-6539EN



**Agilent Technologies**



# The Agilent 5977B GC/MSD and High Efficiency Source (HES) Lowers Detection for Semivolatile Compounds

## Application Note

### Authors

Dale Walker, Melissa Churley, and Harry Prest

### Abstract

The enhanced signal produced by the High Efficiency Source (HES) of the Agilent 5977B GC/MSD allows for flexibility in analytical approach. To assist setting expectations for strategies in semivolatile organic compound analysis, this application note suggests preliminary instrument detection limits for a range of analytes of interest. The results show that picogram or sub-picogram detection in scan mode is possible for a wide range of compounds.



**Agilent Technologies**

## Introduction

Semivolatile organic compounds (SVOCs) are a broad class of environmentally significant contaminants of global interest. These compounds are found on a variety of target analyte lists in GC/MS methods such as the USEPA 8270 and 525 methods, and comparable methods elsewhere. Although listed as targets and appropriate to selected ion monitoring (SIM) in GC/MS analysis, surveying samples by scanning GC/MS provides advantages such as full scan spectra for compound confirmation, tentatively identifying unexpected unknowns in samples that would escape SIM. In the past, scan sensitivity was borderline or insufficient when compared to SIM or the required detection limits. The High Efficiency Source (HES) of the Agilent 5977B GC/MSD represents a revolution in ion source design with greatly enhanced sensitivity that can be exploited to produce scan detection limits for SVOCs that were formerly only approached by SIM. This application note provides preliminary results for instrument detection limits (IDLs) for a few SVOCs across the classes of compounds typical to this analysis.

## Experimental

Many configurations and approaches are used for SVOCs, but common to all is a 5-phase column with a thickness appropriate to the loaded range. Because this was a survey of compounds and a preliminary investigation, an Agilent J&W DB-UI 8270D Ultra Inert GC column (0.25 mm × 30 m, 0.5 μm) was used, as is common in this analysis. This choice would favor higher amounts on-column, although the data suggest that a better approach would be based on a thinner film. Standards were prepared in dichloromethane, and 0.5 μL was injected using a 5-μL syringe in pressure-pulsed splitless mode into a double-taper liner. Replicate injections of a 5 ng/mL standard were used to determine an IDL for each compound. The Agilent 7890B GC with a 5977B HES was operated in scan mode from 50 to 550 u (sampling = 4), and a very low gain factor (0.1) to be pertinent to the desire for a working concentration range. This is essentially a standard configuration, operated in a conservative mode to survey possible IDLs.

## Results and Discussion

Estimated IDLs for the scan data were calculated using the external standard method, and were based on eight consecutive injections of 12 total injections. The average of five IDL determinations was reported as the IDL. As can be seen in Table 1, sub-picogram scan detection is common, with a few compounds showing picogram levels due primarily to lowered compound target ion response. Compound chromatography also played a role in some cases (for example, benzo[b]- and [k]fluoranthene, and so forth).

## Conclusions

Clearly, compound detection in scan mode is now able to discern amounts previously attained only in SIM mode. This advantage allows several analytical strategies to be explored and applied. The shoot less and get more approach means applying split injections with accelerated run times if high concentration levels wish to be maintained. Shooting less sample would also put less matrix in the liner, column, and so forth, and allow the analyst to get more runs before servicing is required. The prep less and save more approach means processing less sample. This would save time and costs not only in collection and transport, but in solvent use and disposal. These dramatically lowered scan IDLs also suggest that SIM IDLs will be enhanced, and so a combination of both strategies is possible to result in the most time and cost effective analysis possible.

Table 1. Agilent 5977B Scan Mode Instrument Detection Limits for SVOCs

<b>Compound</b>	<b>Scan IDL (pg)</b>	<b>Compound</b>	<b>Scan IDL (pg)</b>
Dimethyl phthalate	0.4	<i>o</i> -Cresol	2.5
Diethyl phthalate	1.1	<i>p</i> -Cresol	2.6
Di- <i>n</i> -butyl phthalate	0.9	2,4-Dimethylphenol	0.5
Butyl benzyl phthalate	5.7	2,4-Dichlorophenol	0.5
<i>Bis</i> (2-ethylhexyl) phthalate	0.5	4-Chloro-3-methylphenol	1.6
Di- <i>n</i> -octyl phthalate	2.1	2,4,6-Trichlorophenol	9.1
1,3-Dichlorobenzene	0.3	2,4,5-Trichlorophenol	3.4
1,4-Dichlorobenzene	0.3	Naphthalene	0.2
Benzyl alcohol	3.1	2-Methylnaphthalene	0.4
1,2-Dichlorobenzene	0.3	2-Chloronaphthalene	0.3
1,2,4-Trichlorobenzene	0.3	Acenaphthylene	0.4
Azobenzene	0.6	Acenaphthene	0.8
Hexachlorobenzene	2.1	Dibenzofuran	0.3
<i>Bis</i> (2-chloroethyl) ether	1.4	Fluorene	0.4
<i>Bis</i> (2-chloro-1-methylethyl) ether	0.6	Phenanthrene	0.2
Aniline	0.4	Anthracene	0.3
N Nitroso-di- <i>n</i> -propylamine	2.2	Fluoranthene	0.8
Nitrobenzene	0.4	Pyrene	0.8
4-Chloroaniline	0.9	Benz[a]anthracene	0.4
2-Nitroaniline	1.3	Chrysene	0.3
2,6-Dinitrotoluene	1	Benzo[b]fluoranthene	0.7
3-Nitroaniline	2.8	Benzo[k]fluoranthene	0.7
2,4-Dinitrotoluene	1.4	Benzo[a]pyrene	0.9
4-Nitroaniline	3.8	Indeno[1,2,3-cd]pyrene	0.7
Diphenylamine	0.6	Dibenz[a,h]anthracene	1.3
Phenol	0.60	Benzo[g,h,i]perylene	0.6
2-Chlorophenol	0.5		

## **For More Information**

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)


Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2016  
Printed in the USA  
January 21, 2016  
5991-6612EN



**Agilent Technologies**



# Screen More Drugs with the Agilent GC/MS Toxicology Analyzer with a High Efficiency Source

## Application Note

Forensic Toxicology

### Authors

Melissa Churley and  
Luis Cuadra-Rodriguez  
Agilent Technologies, Inc.

### Introduction

Broad-range screening for drugs in biological samples requires full-spectrum identification confirmation for an unlimited number of targets, as well as spectral identification of nontargets. The Agilent GC/MS Toxicology Analyzer uses Deconvolution Reporting Software (DRS), the Forensic Toxicology Database Library and, when configured with the Agilent 5977B Mass Selective Detector, a high efficiency source (HES). In combination, these technologies screen a greater number of targets at low concentrations while reducing analysis time.



**Agilent Technologies**



## Screen at Lower Concentrations with the High Efficiency Source

Improved screening capability with the HES was demonstrated using human serum, as seen in Table 1. Negative serum, 2 mL, was extracted using the Agilent Bond Elut Certify general drug screen method M2721 [1]. The extract was reconstituted in 0.1 mL methanol and spiked with the GC/MS Toxicology Checkout Mixture (p/n 5190-0471) to yield 10 to 1,000 ng/mL (in vial) of underivatized standards.

Table 1. Lowest injected amount of drugs in spiked serum extract detected with AMDIS using a minimum match factor of 75.

	HES (HES autotune)			Extractor (etune)		
	Minimum amount injected (pg) AMDIS score > 75	AMDIS Score	Equivalent concentration in serum (ng/mL)*	Minimum amount injected (pg) AMDIS score > 75	AMDIS Score	Equivalent concentration in serum (ng/mL)*
Amphetamine	500	94	25	500	75	25
Nicotine	50	92	2.5	50	81	2.5
MDA	500	77	25	500	76	25
MDMA	500	85	25	500	83	25
MDEA	10	76	0.5	500	97	25
Meperidine	10	85	0.5	50	85	2.5
Phencyclidine	50	83	2.5	500	90	25
Methadone	50	87	2.5	500	89	25
Cocaine	50	77	2.5	500	94	25
SKF-525a	50	77	2.5	100	81	5
Codeine	100	88	5	500	90	25
Diazepam	50	90	2.5	50	81	2.5
Hydrocodone	100	91	5	500	90	25
Tetrahydrocannabinol	50	75	2.5	100	78	5
Oxycodone	50	80	2.5	500	83	25
Flunitrazepam	500	88	25	500	75	25
Diacetylmorphine	100	79	5	1,000	83	25
Fentanyl	50	85	2.5	50	77	2.5
Alprazolam	100	76	5	1,000	85	50
Verapamil	50	84	2.5	500	90	25
Strychnine	500	86	25	500	77	25
Trazodone**	> 1,000	(71)	> 50	> 1,000	(68)	> 50

Drugs found at lower concentrations using the HES versus the extractor source are highlighted. Tuning conditions are in parentheses.

\* Assumes 100% recovery from a 2 mL serum sample, reconstitution of extract in 0.1 mL and 1 µL injected.

\*\* The amount of injected trazodone required to achieve a score of 75 exceeds 1,000 pg.

The benzodiazepines oxazepam, lorazepam, temazepam, nitrazepam, and clonazepam were not found at 1,000 pg.

The synthetic opioid fentanyl was added to the checkout mixture. Drug target compounds screened in serum, such as methadone, cocaine, hydrocodone, THC and others, can now be positively identified using full-scan mode at lower concentration (for example, 5 ng/mL for hydrocodone). The HES maximizes the number of ions that are created in the source and transferred into the quadrupole analyzer, which equates to more signal, and thus better sensitivity (Figure 1). This increase in response translates into more drug targets found during the screening process with good library matches.

## NIST Searchable Spectra with Isomer Differentiation

The forensic chemist is required to identify the drug of abuse with the highest degree of scientific certainty. To some extent, the issue of cocaine isomer determination has reinforced the requirement for methods of high specificity [2]. The 5977B GC/MSD with HES allows for the strict spectral integrity that is required to differentiate between cocaine isomers and obtain NIST searchable spectra for detected drugs (Figure 2).




Figure 1. Overlay of TICs of 500 µg standard in serum using A) HES autotune, and B) extractor source etune.




Figure 2. Mass spectrum of 100 µg cocaine spiked in serum (A) compared with NIST spectrum (B). Cocaine is the first hit in NIST. The match factor is 810 (good) [3], and is differentiated from pseudococaine, which has a match factor of 788 (fair). The equivalent concentration is 5 ng/mL based on complete recovery. Excellent NIST library matches ( $\geq 900$ ) for cocaine as a first hit are returned at concentrations above this level.

## Conclusions

The high efficiency source of the Agilent 5977B GC/MSD greatly enhances the signal of drug targets. Resulting spectra are classical and NIST searchable. When combined with Deconvolution Reporting Software, detection levels during screen analysis approach those using SIM mode with derivatization.

## Acknowledgements

The authors thank Bruce Quimby and Fred Feyerherm for productive discussions surrounding this work and also Nathan Contino for his contributions.

## References

1. Anon. *Agilent Bond Elut Certify and Certify II Methods Manual*; Agilent Technologies, Inc. Publication number 5991-4939EN, **2014**.  
<http://www.agilent.com/cs/library/brochures/Bond%20Elut%20Certify%20MethodsManual.pdf>
2. Schlesinger, H. L. Topics in the chemistry of cocaine. *B. Narcotics* **1985**, *1*, 63-78. United Nations Office on Drugs and Crime, Vienna, Austria.  
[https://www.unodc.org/unodc/en/data-and-analysis/bulletin/bulletin\\_1985-01-01\\_1\\_1\\_page006.html](https://www.unodc.org/unodc/en/data-and-analysis/bulletin/bulletin_1985-01-01_1_1_page006.html)
3. Anon. NIST Standard Reference Database 1A, NIST/EPA/NIH Mass Spectral Library (NIST 14) and NIST Mass Spectral Search Program (Version 2.2), User's Guide. National Institute of Standards and Technology, U.S. Department of Commerce, Gaithersburg, MD, USA.  
<http://www.nist.gov/srd/upload/NIST1aVer22Man.pdf>

## For More Information

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

For Forensic Use only.  
Information is subject to change without notice.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

© Agilent Technologies, Inc., 2015  
Printed in the USA  
October 14, 2015  
5991-6294EN



**Agilent Technologies**



The Agilent 5977B GC/MSD

# THE SOURCE OF NEW POSSIBILITIES

The Measure of Confidence



**Agilent Technologies**

# HIGH EFFICIENCY, HIGH CAPABILITY

## Agilent's new High Efficiency Source (HES) gives you 10x improvement in capabilities.


The Agilent 5977B HES GC/MSD System builds on a long tradition of providing the most trusted single-quadrupole GC/MS system. We've now pushed the technological frontier even further by introducing our newest advancement: the HES – High Efficiency Ion Source.

Fundamentally, mass spectrometry sensitivity depends on the number of ions measured. The ultra-efficient EI source of the 5977B HES GC/MSD maximizes the number of ions that are created and transferred out of the source body and into the quadrupole analyzer. This novel design revolutionizes single-quadrupole MS performance, offering two distinct advantages:

### Advantage #1

#### 10x Greater Sensitivity

Take advantage of the improved sensitivity and see detection limits as low as 1.5 fg IDL – performance that brings yesterday's triple-quadrupole limits of detection into today's single-quadrupole laboratory.



### Advantage #2

#### 10x Less Sample Required

You'll spend less time performing sample prep and maintenance while saving on shipping costs at the same time.



#### Analyzers for guaranteed chromatographic performance

The exceptional performance of the GC/MSD systems is validated by data from common food safety, environmental, and forensics toxicology methods.



#### Go green, go dry with the IDP-3 oil free Vacuum Pump

It offers a quieter laboratory environment, no oil contamination and lower cost of ownership compared to standard oil rotary vane pump.

# OPTIMIZE WORKFLOWS IN FUNDAMENTALLY NEW WAYS

## Work smarter with integrated GC, MS, and software technologies

Agilent has completely reimagined electron ionization, allowing you to reimagine the way you think about your laboratory's operational workflow.

See analytes that have been hidden up to now. Inject a fraction of the valuable sample that you currently require. Collect, ship, store, and process a fraction of the per-sample volume that you have in the past. Handle and dispose of less hazardous solvent. The list of potential operational efficiencies is endless!



## Agilent CrossLab

From Insight to Outcome

End to end laboratory support from the Agilent service experts to keep your instruments running at peak performance with minimal downtime. We offer solutions that include repair, preventative maintenance, compliance verification, education, and a host of other services tailored to meet your needs. Ask us how we can support your laboratory today!



### Smart technology aligns GC and MS operation

The Agilent 7890B GC – with its efficient protocols and fully synchronized MS operation – is a dynamic partner for Agilent GC/MSD systems.

### Integrated software tools simplify method development

From instrument settings to data analysis and reporting, MassHunter puts you in control – and makes your GC/MSD analysis routine.

As available sample becomes smaller, increasingly active, and more complex, you simply cannot afford interferences caused by flow path activity.

**Agilent Inert Flow Path solutions** ensure a reliably inert flow path for higher sensitivity, accuracy, and reproducibility.



Learn more

**[www.agilent.com/chem/5977B](http://www.agilent.com/chem/5977B)**

GC column selection tool

**[selectgc.chem.agilent.com](http://selectgc.chem.agilent.com)**

Inert flow path

**[www.agilent.com/en-us/  
promotions/inertflowpath](http://www.agilent.com/en-us/promotions/inertflowpath)**

Buy online

**[www.agilent.com/chem/store](http://www.agilent.com/chem/store)**

Find a local Agilent customer center  
in your country

**[www.agilent.com/chem/contactus](http://www.agilent.com/chem/contactus)**

USA and Canada

**1-800-227-9770**

**[agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)**

Europe

**[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)**

Asia Pacific

**[inquiry\\_lsca@agilent.com](mailto:inquiry_lsca@agilent.com)**

For research use only. Information, descriptions,  
and specifications in this publication are subject  
to change without notice. Agilent Technologies  
shall not be liable for errors contained herein  
or for incidental or consequential damages in  
connection with the furnishing, performance,  
or use of this material.

© Agilent Technologies, Inc., 2015  
Printed in the USA, October 15, 2015  
5991-6299EN



**Agilent Technologies**